



## MyCare<sup>TM</sup> INSITE

### Quantitative Antipsychotic Blood Levels at the Point of Care

#### Designed to be **Where** and **When** it's **Needed Most**

- Allows near patient testing
- Rapid actionable results (< 10 minutes)
- Small, compact design
- Simple, intuitive operation
- CE Marked and Health Canada approved
- Capillary blood sampling
  - Patient-friendly
  - Collection can be completed by trusted care team member
  - Requires no phlebotomist
  - Safer for staff and patients

#### MyCare Insite Test Menu

Clozapine\*  
Total Aripiprazole\*  
Total Risperidone  
Quetiapine  
Olanzapine

\*CE Marked

Learn more at  
[MyCareTests.com](http://MyCareTests.com)



# Simple, Accessible, and Affordable

MyCare<sup>™</sup>  
INSITE

The MyCare<sup>™</sup> Insite improves access to accurate diagnostic testing by providing antipsychotic blood levels at the point of care. Healthcare providers can now make rapid triage and treatment decisions when diagnosing a patient's condition or monitoring treatment response. Compared to sending samples to a laboratory, the MyCare Insite provides numerous clinical and operational advantages while also providing significant cost savings. The Insite's simplicity of operation, its compact design and rapid actionable results improves the accessibility and clinical benefit of antipsychotic blood levels.

- Maintain therapeutic blood levels to avoid potential mental status deterioration and minimize adverse effects
- Enable timely responses to changes in clinical status
- Reduce the length of stay in specialized environments
- Improve patient experience and engagement with immediate test results and feedback

## 4-Step Procedure

### 1 Collect patient sample




### 2 Insert the cartridge



### 3 Run the sample



### 4 View patient result



## MyCare Insite Features

### Compact, Portable Design

Dimensions: 16 x 13 x 14.5 cm

Weight: 2.2kg

### Easy to Operate

Intuitive software similar to mobile devices

Software available in 17 languages

Maintenance-free

Factory calibration

### Connectivity

The MyCare Insite operates in three modes:

- Standalone
- EHR connected (LIS or HIS compatible)
- Optional cloud-based services provide capabilities to remotely view test results

